


# A VERSATILE COMBINATION

New products require a highly flexible approach in fill and finish. Schott, ARaymond and groninger offer solutions for future challenges in drug manufacturing and provide maximum flexibility to customers.



## SCHOTT AG

Gregor Deutsche  
Director Business  
Development Sterile Solution  
Phone +49 6131 66 - 4049  
gregor.deutsche@schott.com

## ARaymondlife

Pascal Sircoulomb  
Business Development Director  
Phone +33 456 52 - 6175  
pascal.sircoulomb@araymond.com

## groninger & co. gmbh

Marius Römer  
Product Manager  
Business Development  
Phone +49 7951 495 - 3259  
m.roemer@groninger.de

# FLEXIBILITY

for your next batches



## FLEXIBLE IN PRIMARY PACKAGING


### Schott iQ®



Develop, deliver and improve injectable drugs like never before with the SCHOTT iQ® holistic platform to standardize ready-to-use (RTU) packaging components and improve the quality of RTU containers.

## FLEXIBLE IN CLOSURE

### RayDyLyo®




RayDyLyo® plastic push fit cap, a new standard for vial capping, that replaces aluminum caps



Suitable for ISO 8362-2 rubber stoppers, pre-validated by the customer

Suitable for ISO 8362-1 vials


Packaged in standard bag and RTP bag:



## FLEXIBLE IN PROCESSING

### Filling and closing machine

Output range: 4,800 – 60,000 objects / h



groninger  
PHARMA

### Highlights:

- Process security using closing force measurement designed for push fit caps
- Flexibility in processing of all nested components just by changing format parts
- Scalability by applying the same processes for each output range